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en

1. Product Description/Application Area

Splint-Fix is a fixing and positioning splint for the immobilization of following joints: finger joints, wrists, elbows and ankles.

2. Indications

The fixing and positioning splint is applied for the support of e.g. smooth peripheral infusion therapies and/or invasive blood pressure monitorings.

3. Contraindications

Splint-Fix **MAY NOT** be used for fractures, luxations, open wounds/skin defects and circulatory disturbances at the place of its application.

4. Patients Targeted

For patients with peripheral venous and arterial accesses where a smooth function is only provided by immobilization.

5. Intended Users

Nursing staff, medical professionals, nursing relatives

6. Application Description

1 | Form the splint beforehand in the position desired and place it, depending on the indication, to the interior or exterior of the joint to be immobilized.

2 | Tie the elastic straps around the joint respectively particular extremity.

3 | Close the straps with the hook and loop fastener; to improve stabilisation also crosswise.

7. LOT Number

The LOT number is imprinted on the product and outer packaging.

8. UDI Carrier

The UDI-Code is imprinted on the product and outer packaging.

9. Expiry Date

The product is not subject to any expiry date.

10. Storage Conditions

Keep Splint-Fix dry and away from sunlight.

11. Handling Conditions

In case of damaged packaging or inappropriate storage the product **MAY NOT** be used and has to be destroyed.

12. Sterile Condition

The product is **NON-STERILE**.

13. Sterilisation Process

Not applicable

14. Sterile Packaging

Not applicable

15. Warning Notices

Application during MR examination

With consideration of the following conditions the application in MR is **partly** possible:

MR conditional ▲

Non-clinical testing has demonstrated that the Splint-Fix fixing and positioning splints are MR conditional. A patient with this device can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 12.900 G/cm (129 T/m)
- **Maximum force product** of 236.000.000 G²/cm (236 T²/m)
- **Theoretically estimated** maximum whole body averaged (WBA) specific absorpti-on rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, Splint-Fix 22 is expected to produce a maximum temperature rise of less than 3,6 °C (2 W/kg, 1,5 Tesla) **RF-related** temperature increase with a background temperature increase of ≈ 0,7 °C (2 W/kg, 1,5 Tesla) and Splint-Fix 13 a maximum temperature rise of less than 1,0 °C (2 W/kg, 3 Tesla) **RF-related temperature increase with a background temperature increase** of ≈ 1,0 °C (2 W/kg, 3 Tesla) after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 2,3 mm from Splint-Fix 9 when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

ATTENTION: MR examinations with different settings can cause serious health damages in patients. The elastic straps **MAY NOT** be applied too tight (risk of circulatory disturbances).

16. Single Use Notice

The product is a disposable product.

17. Notice on Reuse

DO NOT reuse. The reuse of disposable products holds a potential risk for patients and users. This can result in contamination and/or effects on the functioning. Contamination and/or limited functioning can cause injuries, diseases and the death of a patient.

18. Reprocessing Notice

The product **MAY NOT** be reprocessed.

19. Notice on Processing Methods

Not applicable

20. Disposal

The product can be disposed with normal clinical waste or according to local regulations. Special guidelines of medical facilities for the disposal (of e.g. multi-resistant pathogens MRSA, VRE) need to be taken into consideration.

21. Notice on the Application with other Products

Not applicable

22. Notice on Medical Devices/Custom-made Products/Products for Clinical Testings

Medical device of Class I.

23. Notice on the Application by Laymen

Splint-Fix is easy to apply and does not require special abilities or equipment. In case of unexpected complications immediately inform a doctor.

24. Notice on Reporting Requirements

Please inform the manufacturer and the responsible authority of the member state in all cases of serious incidents in connection with the product.

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